ENGLISH

CAUTION: Refer to instruction manual/booklet(IFU).

NOTICE TO USER: In the event a serious incident has occurred in relation to the use of this device and/or a device deficiency occurs, please notify your local representative and/or manufacturer as soon as possible via email to complaints@angiosafe.com. Notify the applicable competent authority in which the user is established. Refer to the Return of Devices section of the IFU.

Contents

One (1) Santreva™-ATK Endovascular Revascularization Catheter, REF: AS-1963-02

Intended Purpose/Use

The Santreva-ATK endovascular revascularization catheter is intended for crossing stenotic lesions, including Chronic Total Occlusions (CTOs), together with plaque compression and channel formation in peripheral arteries (e.g., femoropopliteal arteries) to facilitate placement of a guidewire in the true lumen after crossing and prepare the vessel for further treatment.

Indications for Use

Santreva-ATK is an endovascular revascularization catheter intended to facilitate placement of a guidewire in the true lumen of peripheral vessels (e.g., femoropopliteal arteries) after crossing CTOs in patients with Peripheral Arterial Disease (PAD) and to prepare the peripheral vessel for further treatment with other interventional devices per the physician's discretion.

Intended Users and Environment

This product is designed for use by physicians trained in and familiar with percutaneous interventional techniques in a fully equipped catheterization laboratory or vascular surgery suite.

Intended Patient Population

Adults, 22 years of age or older. To reduce the risk of vessel damage, the Santreva-ATK endovascular revascularization catheter should be used in subjects with target lesion reference vessel diameter(s) of ≥ 3.0 mm and ≤ 10 mm by visual estimate.

Device Description

The Santreva-ATK (Figure 1, Figure 2, and Figure 3) is an endovascular revascularization catheter, which is intended for intraluminal crossing, together with plaque compression and recanalization of chronic total occlusions (CTOs) in peripheral vessels (e.g., femoropopliteal arteries) to prepare the vessels for final treatment. The catheter consists of an inner single lumen stainless steel torque cable concentric to an outer braided shaft with polymer jacket. The flexural and axial properties of the Santreva-ATK endovascular revascularization catheter have been designed for femoral arterial vessel access and delivery to the CTO disease location. The inner lumen of the catheter torque cable is designed to facilitate placement of a guidewire into the true vessel lumen after CTO crossing and recanalization.

The distal end of the catheter has a tapered stainless steel cutting tip that measures 1.35 mm (0.053 inch) in outer diameter. This cutting tip has an integrated 0.20 mm (0.008 inch) outer diameter stainless steel wire cutting loop mounted at the distal face to puncture, displace, and compress the plaque as the mechanism of CTO crossing.

The cutting tip and loop combination is maintained within the lumen of the artery and CTO body during crossing by an electropolished nickel titanium (nitinol shape set memory metal) centering system component of 2.4 mm (0.095 inch) outer diameter and 7.6 mm (0.300 inch) length. Together, the cutting tip and centering system have a combined length of 9.8 mm (0.386 inch). The centering system has three shape set wings that are laterally flat when expanded which glide safely in the open lumen of the artery and within the CTO body during crossing. At their attachment junctions to the distal catheter shaft, the collapsible centering system wings are sloped in each direction to present a tapered interface to the CTO body plaque and with accessories such as guide catheters or guiding sheaths during both advancement and retraction of the catheter.

The Santreva-ATK endovascular revascularization catheter user interface consists of a palm-sized handle with an integrated stainless steel and manually controlled handle wheel connected to the inner torque cable to provide cutting tip and loop rotation. During use, the operator grasps the outer catheter shaft to advance the device in the distal direction while rotating the handle wheel such that the cutting tip and loop combination punctures the entrance cap to the CTO body plaque. As the operator continues to advance the catheter shaft and manually rotate the handle wheel, the connected torque cable translates the rotation to the cutting tip and loop ensemble such that the CTO plaque is simultaneously punctured by the rotating cutting loop, radially displaced and compressed by the tapered slope of the rotating cutting tip. As the user advances the device through the CTO body, the plaque is further compressed laterally by the centering system wings, thereby crossing and recanalizing the CTO intraluminally in a controlled manner. Once the cutting tip and loop of the catheter is advanced to the exit cap of the CTO body, the user advances a guidewire tip through the channel formed and into the distal lumen beyond the CTO and the Santreva-ATK catheter is retracted from the peripheral vasculature while leaving the guidewire in place for use in final treatment.

The guidewire is placed using the inner single lumen of the Santreva-ATK catheter which extends from an integrated luer fitting at the proximal end of the handle through the torque cable to a 0.43mm (0.017 inch) inner diameter exit port of the cutting tip. The catheter is compatible with 0.36 mm (0.014 inch) diameter guidewires.

The Santreva-ATK endovascular revascularization catheter is compatible with 6 Fr (0.070 inch / 1.78 mm Minimum ID) or larger guide catheters and 5 Fr (0.070 inch / 1.78 mm Minimum ID) or larger guiding sheaths when not using a guide catheter. The catheter working length is 135 cm from the handle strain relief to the distal tip of the device.

Figure 1: Santreva™-ATK Endovascular Revascularization Catheter Handle, Catheter Shaft, Expanded Centering System and Distal Tip


Figure 2: Santreva™-ATK Handle and Catheter Shaft Lengths



Figure 3: Santreva™-ATK Expanded Centering System and Distal Tip

Contraindications

- The device is not intended for use in the renal, mesenteric, cerebral, or coronary vasculature.
- Contraindicated for use in patients with known hypersensitivity to nickel or titanium.
- Contraindicated for use in patients who cannot tolerate anticoagulant or anti-platelet therapy.
- Contraindicated for use in patients with unresolved bleeding disorders.

Warnings

- Do not use the Santreva-ATK catheter with hydrophilic coated guidewires.
- Do not use in a target lesion located in a bypass graft or stent (i.e., in-stent restenosis).
- Before the insertion of the device, administer appropriate anti-coagulant.
- The device is intended for single patient use only. Do not re-sterilize or reuse, as this can result in compromised device performance.
- Do not use if the package is opened or damaged.
- Do not use after the “Use By” date on the product label.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or government policy.

Precautions

- Do not use without completely reading and understanding this document.
- Store device at ambient conditions in a catheterization lab or storage room.
- Do not expose the catheter to organic solvents (e.g., alcohol).
- The outer box is not a sterile barrier. The pouch contained within the outer box is the sterile barrier. Only the contents of the inner pouch are sterile. The outside surface of the inner pouch is not sterile. Do not remove the contents of pouch until immediately prior to the use of device.
- Excessive bending or kinking of the catheter may affect performance.
- Torquing the catheter excessively may cause damage to the product.
- If the catheter is believed to become kinked, unintentionally damaged, and/or structural integrity of the device is thought to be impacted, withdraw the catheter and replace with a new unit to continue with the procedure.
- If strong resistance is felt during manipulation, determine the cause of the resistance before proceeding further. If the cause cannot be determined, withdraw the catheter.
- If the device is withdrawn during the CTO crossing procedure, inspect the device tip and catheter shaft for damage prior to reinsertion into the patient.

Expected Clinical Benefits

The expected clinical benefits of Santreva-ATK when used as intended in peripheral vessels (e.g., femoropopliteal arteries) of patients with Peripheral Arterial Disease (PAD) include:

1. Crossing stenotic lesions, including Chronic Total Occlusions (CTOs).
2. Plaque disruption and compression without risk of embolism, flow-limiting dissection, and major perforation.
3. Formation of an intraplaque and angiographically visible channel.
4. Facilitation of guidewire placement in distal true lumen through the device.
5. Leaving the lesion and vessel prepared one step for further imaging and final treatment.

One-step vessel preparation with Santreva-ATK enables further endovascular treatment and can also lead to shorter procedure times. These potential clinical benefits of the Santreva-ATK device are intended to be achieved with a high degree of technical success and low occurrence of periprocedural complications.

Summary of Clinical Safety and Performance

The RESTOR-1 study titled “Safety and Effectiveness Study of the AngioSafe Santreva™-ATK Endovascular Revascularization Catheter (RESTOR-1),” was a prospective, single-arm, multi-center, pivotal clinical study designed to enroll subjects with a prior diagnosis of Peripheral Arterial Disease (PAD) and peripheral Chronic Total Occlusion (CTO) in the peripheral arteries of lower limbs.

The RESTOR-1 study enrolled a total of 179 subjects in the pre-screening portion of the study at a total of 14 investigational sites, all of which were located in the United States (U.S.). A total of 132 subjects were eligible for the Baseline Visit Screen, and of these, 91 subjects were eligible for the study procedure. A sample size of 79 subjects had post-procedure data reviewed by the Core Laboratory (Full Analysis Set, FAS); this 79-subject population was used for safety analyses (Initial FAS and Safety Set n=79). Five (5) subjects were removed from the FAS due to protocol changes resulting in a final FAS population of 74 subjects (FAS n=74) used for effectiveness analyses. Two (2) additional subjects were removed from the FAS for major protocol violations resulting in a Per Protocol (PP) analysis set of 72 subjects (PP n=72).

Study Endpoints	
Primary Effectiveness Endpoint	Clinical success, defined as the ability to facilitate placement of a guidewire into the distal true lumen of a femoropopliteal artery with CTO in the absence of device-related MAEs through discharge or 24-hours post-procedure, whichever occurred sooner.
Secondary Endpoints	<ol style="list-style-type: none"> 1. Technical success, defined as the ability of the study device to facilitate placement of a guidewire into the distal lumen. 2. Procedural success, defined as Technical Success without a procedural complication within 30 days after the procedure. 3. Evaluation of intraluminal CTO crossing (as assessed by intravascular ultrasound, IVUS). 4. Primary endpoint in the subgroup stratified by the degree of calcification.
Safety Endpoints	<ol style="list-style-type: none"> 1. Frequency of device-related MAEs through discharge or 24-hours following the procedure, whichever occurred sooner. 2. Frequency of MAEs through Day-30. MAEs were defined as: all-cause death, adverse events leading to unplanned amputation of the treated extremity or unplanned endovascular and/or surgical revascularization of the treated extremity, and distal embolization requiring additional treatment after crossing the lesion with the Santreva-ATK device.

Primary inclusion criteria were peripheral artery disease (PAD), with Rutherford Clinical Classification 2 – 5 (“moderate claudication” to “minor tissue loss”) with a Chronic Total Occlusion (CTO) in the peripheral arteries of the lower limbs. Confirmation of stenosis was completed by angiography or Duplex ultrasound. Subjects who presented with Acute Limb Ischemia (ALI), previous major amputation above the ankle in the extremity that was to be treated were excluded from the study. Subjects were assessed at pre-treatment, treatment, and hospital discharge. Follow-up visits were scheduled for Day-30. The RESTOR-1 study utilized an independent angiographic Core Laboratory, as well as a Clinical Events Committee (CEC), to assess the outcomes data and adjudicate adverse events.

The null hypothesis was that the AngioSafe Santreva-ATK Endovascular Revascularization Catheter would perform better than 0.70 - the lower bound of the 95% Confidence Interval (CI) of the technical success rate. The primary effectiveness endpoint of this study, the clinical success of the Santreva-ATK Endovascular Revascularization Catheter, was met and was achieved in 87.8% of the study participants (65 subjects) in the FAS population, with a 0.78 lower bound of the 95% CI, with over 70% of the subjects having CTOs with moderate to severely calcified plaques. In the PP population, the success rate was 88.9%, with a 0.79 lower bound of the 95% CI with the same level (over 70%) of moderate to severely calcified plaques. Note that in the PP population, based on further analysis, the Santreva-ATK device was delivered to non-CTO vessels in two subjects. In this case, Santreva-ATK successfully facilitated wire placement in the distal true lumen in 63 of 70 subjects, achieving a 90% final success rate, with over 70% of the subjects having CTOs with moderately to severely calcified plaques.

The secondary endpoints of technical and procedural success were achieved at the same rate as the primary endpoint of clinical success since procedural success was defined in this study as technical success without a procedural complication within 30 days after the procedure and there were no procedural complications within 30 days after any of the procedures. Clinical success of the primary endpoint was stratified by levels of calcification, which included two groups: none to mild calcification and moderate to severe calcification. Results are consistent across all levels of calcification. The success rate for none to mild calcification was 90.4%, and for the moderate to severe calcification 86.7%. Average treated CTO length was 131.6 mm, average target lesion diameter was 5.7 mm, target lesion severe calcification mean was 34.1%, median crossing time was approximately 9 minutes, and mean crossing time was approximately 25 minutes.

The AngioSafe Santreva-ATK Endovascular Revascularization Catheter created an angiographically visible, and IVUS-confirmed, track in all RESTOR-1 subjects (100%) who had technical success and for whom angiographic and IVUS data was collected after CTO crossing completion. The average track formed by the AngioSafe Santreva-ATK device was approximately 2.8 mm in diameter, which was measured by quantitative angiography or quantitative IVUS. The centering system comprises the largest feature of the Santreva-ATK device, with a diameter of 2.4 mm. Therefore, an average of approximately 2.8 mm track diameter formed through the traversed CTO is reflective of Santreva-ATK's device course through the lesion. This newly formed intraluminal track resulted in an average lumen gain (calculated by dividing the mean track diameter with mean reference vessel diameter) of 59%.

There were no device-related Major Adverse Events (MAEs) through discharge or 24-hours following the procedure, whichever was sooner. There was one (1) MAE within the 30 days following the procedure, which was neither procedure-related nor device-related.

Potential Complications

The use of this product carries the risks associated with peripheral vascular angioplasty, including thrombosis, vascular complications, and/or bleeding events. The risks associated with standard PTA procedures are reported in the published literature and include the following:

<ul style="list-style-type: none"> • Access site pain • Allergic reaction to contrast medium, anticoagulant, antithrombotic therapy, or device materials • Aneurysm • Arrhythmias • Arterial dissection • Arterial perforation • Arterial rupture • Arterial spasm • Arteriovenous fistula • Bleeding complications • Concomitant medication complications (drug reactions, bleeding from antiplatelet/anti-coagulation agents) • Death • Emboli (air, tissue, thrombus or atherosclerotic emboli) • Emergency or non-emergency arterial bypass surgery • Entry site complications • Fever 	<ul style="list-style-type: none"> • Fracture of the guidewire or any component of the device that may or may not lead to device embolism, serious injury or surgical intervention • Hematoma • Hemorrhage at the vascular access site • Hemolysis • Hypertension • Hypotension • Infection • Ischemia • Myocardial infarction • Pseudoaneurysm • Renal failure • Restenosis of the treated segment • Sepsis • Shock/pulmonary edema • Thrombosis • Total occlusion of the peripheral artery • Vascular complications which may require surgical repair (conversion to open surgery)
--	---

Guidewire and Support Accessory Compatibility

- 0.014 inch / 0.36 mm Nominal OD Guidewire (Note: Do not use the Santreva-ATK catheter with hydrophilic coated guidewires)
- Minimum size of 5 Fr (0.070 inch / 1.78 mm Minimum ID) or greater Introducer or Guiding Sheath **without a Guide Catheter**.
- Minimum size of 6 Fr (0.070 inch / 1.78 mm Minimum ID) or greater Guide Catheter **with 6 Fr or greater Introducer or Guiding Sheath**

Note: Use of a 6 Fr or greater size Guide Catheter is recommended for maximum support.

Device Preparation

1. Use sterile technique to carefully remove the device from the packaging. Inspect to ensure the device exhibits no signs of damage.
2. Flush the Santreva-ATK catheter lumen with heparinized saline using the proximal luer of the handle.
3. If desired, pre-load a 0.014" guidewire into the proximal insertion luer of the handle and advance the guidewire until it is approximately 1 cm proximal to the distal tip of the device.
Note: The guidewire may be backloaded into the device tip as necessary for over-the-wire technique during the procedure.
4. If desired, pre-load the Santreva-ATK catheter into the support accessory (selected sheath or guide catheter) prior to inserting into the vasculature by advancing the device tip until just proximal of the support accessory tip.
Note: Rotate the Santreva-ATK catheter tip and centering system counterclockwise during insertion into the proximal hub of the support accessory to assist the folding and compression of the centering system.

Directions for Use

1. Prepare the insertion site using sterile technique.
2. Achieve vascular access and place an introducer or guiding sheath.
3. Inject a weight-based Unfractionated Heparin (UFH) bolus dose followed by as needed additional boluses to achieve and maintain an Activated Clotting Time (ACT) of 250-400 sec.
4. Insert the Santreva-ATK catheter into the selected support accessory or a pre-loaded assembly of Santreva-ATK catheter/guidewire/support accessory into the selected introducer and/or guiding sheath.
5. Advance the selected support accessory tip as close to the entry cap of the CTO for maximum support.
6. If not pre-loaded, insert a 0.014" / 0.36 mm guidewire into the proximal insertion luer of the handle and advance the guidewire until it is approximately 1 cm proximal of the distal tip of the Santreva-ATK catheter.
7. Advance the Santreva-ATK catheter tip out of the support accessory tip until the nitinol centering system has exited and expands.
8. With the guidewire remaining within the Santreva-ATK catheter proximal of the distal tip, advance the Santreva-ATK catheter tip to the entry cap of the CTO.
9. Advance the Santreva-ATK catheter shaft with the guidewire remaining in place with one hand to allow the distal tip to penetrate the entry cap of the CTO while simultaneously rotating the torque wheel on the handle with the other hand in either a clockwise or counterclockwise (or alternate directions) such that the cutting loop at the tip of the device facilitates cap penetration.
10. Advance the Santreva-ATK catheter tip and centering system with the guidewire remaining in the Santreva-ATK catheter through the CTO by advancing the catheter shaft with one hand while continuing to rotate the torque wheel on the handle with the other hand in either a clockwise or counterclockwise direction (or alternate directions) until the distal tip is approximately 1 cm from the exit cap of the CTO.
11. Cautiously continue to advance the Santreva-ATK catheter and distal tip while rotating torque wheel up to the distal cap until the tip and cutting loop with the guidewire remaining in the Santreva-ATK catheter are through the exit cap.
Note: Only the Santreva-ATK distal tip and cutting loop should cross the exit cap to minimize embolization risk.
12. Advance the 0.014" / 0.36 mm guidewire distally out of the Santreva-ATK catheter tip past the exit cap of the CTO into the distal vessel lumen.
13. While stabilizing the guidewire distal to the CTO, carefully withdraw the Santreva-ATK catheter from the CTO using over-the-wire technique, while confirming that the guidewire tip remains in the true vessel lumen beyond the CTO.
14. Retract the Santreva-ATK catheter tip and centering system into the tip of the support accessory using simultaneous axial pull and counterclockwise rotation of the Santreva-ATK catheter shaft.
15. Continue retracting the Santreva-ATK catheter from the support accessory to remove the catheter from the patient and proceed with prescribed treatment, if required.

Return of Devices

1. If any portion of the AngioSafe Santreva-ATK catheter fails, is unintentionally damaged, and/or structural integrity of device is thought to be impacted prior to or during a procedure, immediately discontinue use of this device and notify the device manufacturer via email at complaints@angiosafe.com.
2. If any Adverse Event occurs in relation to the use of this device, notify the manufacturer via email at complaints@angiosafe.com and your local representative.
3. In the event a serious incident has occurred in relation to the use of this device, notify the applicable regulatory authority of the country in which the user is established.

Symbol Key

Symbol	Description	Symbol	Description
	Catalog number		Lot / Batch number
	Date of Manufacture		Use-by Date
	Packaging Unit		Medical Device
	Refer to instruction manual/booklet (IFU)		Do not use if packaging is damaged
	Caution		Non-Pyrogenic
	Keep dry		Do not re-use
	Do not resterilize		Single sterile barrier system with protective packaging outside
	Sterilized using radiation		Manufacturer
	Unique device identifier		Authorized representative in the European Community/European Union
	Authorized representative for Switzerland		Authorized representative in the UK
	Product meets the European Union's (EU) health, safety, and environmental protection standards		

 <p>AngioSafe, Inc. 5215 Hellyer Ave, Suite 240 San Jose, CA 95138, USA TEL: (408) 217-5480</p> 	 OMC NI MEDICAL LIMITED Unit - 763, Moat House Business Centre 54 Bloomfield Avenue Belfast, Northern Ireland, BT5 5AD
 OMC Medical Limited Planet House, North Heath Lane Horsham, West Sussex, RH12 5QE, UK	 OMCS Medical GmbH Aegeristrasse 5 Zug 6300, Switzerland

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